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Reasons for revision	n	per 100 OCY
Infection	21	0.27
Mechanical loosening	18	0.23
Pain	14	0.18
Fracture of bone around the joint	8	0.10
Implant instability	6	0.08
Breakage of the implant	2	0.03
Wear	0	0.00
Metal effects	0	0.00
Other	16	0.21

1. Scott CE, et al. *J Bone Joint Surg Br.* 2010; 92-B(9):1253-1258.
2. Bourne RB, et al. *Clin Orthop Relat Res.* 2010; 468:57-63.
3. Huch K, et al. *Ann Rheum Dis.* 2005; 64:1715-1720.

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Linking Kinematic Theories to Implant Design

In order to understand the functional outcome differences produced by various implant designs, it is critical to first examine the thought progressions that led to these designs. Many of today's popular designs still adhere to outdated kinematic theories that require less conforming tibiofemoral articulation to produce rollback. However, instead of producing rollback, these designs have been shown to produce instability and limited function.^{2013,25} The eMP™ Knee System restores function by replicating the natural AP stability and medial-pivot motion of the knee.²⁰¹⁶

Kinematic Theories

First Generation

"Four-Bar Link" Theory:

- ▶ "The kinematics of the knee is guided by a four-bar link."⁹
- ▶ Axis of rotation located at the crossing point of the cruciate ligaments.
- ▶ Demonstrated as progressive rollback of the femur relative to the tibia.

1970

"Four-Bar Link Theory"

Implant Designs

Multi-Radius Knees:

- ▶ Designed to produce rollback.
- ▶ Rollback dictates limited tibiofemoral conformity.
- ▶ Commonly associated with anterior sliding, instability, and anterior knee pain.

Kinematic Theories

Second Generation

Single-Axis of Rotation:

- ▶ Hollister AM. The axes of rotation of the knee. 1993.¹
- ▶ Single flexion – extension axis runs through collateral ligament origins.
- ▶ Longitudinal axis of rotation passes through insertion of ACL on tibial plateau.

1990

Single-Axis of Rotation

Implant Designs

Single-Radius Knees:

- ▶ Designed to produce symmetric rollback and freedom of rotation about longitudinal axis.
- ▶ Lack of tibiofemoral congruency to allow central rotation.

Kinematic Theories

Latest Generation

Medial-Pivoting Kinematics:

- ▶ Medial side of the knee stays stable and demonstrates spinning kinematics.
- ▶ Lateral side moves anterior and posterior to accommodate rotation of the tibia.
- ▶ Sagittal view shows circular articulating surfaces.

2000

Medial-Pivoting Kinematics

Implant Designs

The Medial-Pivot Knee:

- ▶ Designed to promote medial stability and lateral mobility.
- ▶ Constant radius femoral component, both in the sagittal and coronal plane, is matched with a highly conforming medial compartment and less conforming lateral compartment to allow translation.

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
An anatomical illustration of a human lower body, showing the legs and feet. The knees are replaced with a mechanical system called eMP. The background is a blue gradient with a large, faint, stylized eye shape. The text 'eMP' is written in a large, white, serif font in the upper left corner.

eMP™

Evolution®
Medial-Pivot Knee System
The ACL-PCL Substituting Knee™

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Unicompartmental Knee Arthroplasty Provides Significantly Greater Improvement in Function than Total Knee Arthroplasty Despite Equivalent Satisfaction for Isolated Medial Compartment Osteoarthritis

Casper

J Arthroplasty. 2019 Jun 26. pii: S0883-5403(19)30634-5. doi: 10.1016/j.arth.2019.06.043. [Epub ahead of print]

The Combination of Inlay Patellofemoral Arthroplasty and Medial Unicompartmental Knee Arthroplasty Versus Total Knee Arthroplasty for Mediopatellofemoral Osteoarthritis: A Comparison of Mid-Term Outcomes

Uluyardimci E¹, Isik C², T

Acta Orthop Belg. 2018 Sep;84(3):237-244.

Bi-unicompartmental versus total knee arthroplasty: long term results.

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J Arthroplasty. 2014 Jan;29(1):48-51. doi: 10.1016/j.arth.2013.04.024. Epub 2013 May 18.

Joint awareness in different types of knee arthroplasty evaluated with the Forgotten Joint score

Thienpont E¹, C

Knee Surg Sports Traumatol Arthrosc. 2017 Mar;25(3):681-686. doi: 10.1007/s00167-015-3868-1. Epub 2015 Nov 21.

Unicompartmental knee arthroplasty versus total knee arthroplasty: Which type of artificial joint do patients forget?

Zuiderbaan H J Arthroplasty. 2017 May;32(5):1453-1459. doi: 10.1016/j.arth.2016.11.034. Epub 2016 Nov 27.

Differences in Patient-Reported Outcomes Between Unicompartmental and Total Knee Arthroplasties: A Propensity Score-Matched Analysis.


Kim M

RESULTS: One-hundred and thirty patients were included. Sixty-five patients underwent medial UKA and 65 patients underwent TKA. At both follow-up points, the FJS was significantly higher in the UKA group (FJS 1 year 73.9 ± 22.8 , FJS 2 year 74.3 ± 24.8) in contrast to the TKA group (FJS 1 year 59.3 ± 29.5 ($p = 0.002$), FJS 2 year 59.8 ± 31.5 , ($p = 0.004$)). No significant improvement in the FJS was observed between 1- and 2-year follow-up of the two cohorts.

CONCLUSION: Patients who undergo UKA are more likely to forget their artificial joint in daily life and consequently may be more satisfied.

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




Design and kinematics in total knee arthroplasty.

Digennaro V¹, Zambianchi F, Marcovigi A, Mugnai R, Fiacchi F, Catani F.

CONCLUSIONS: Journey BCS showed statistically significant better KOOS results. The higher posterior femoral rollback observed in the kinematic assessment of this design, associated with a better patellofemoral design, may be the reason for better clinical outcome. The reported cases of stiffness and anterolateral joint pain could be attributed to excessive medial and lateral tibiofemoral posterior translation. The NRG group demonstrated good axial rotation, but this was not coupled with physiological kinematic patterns. Patellofemoral pain can be explained by a less friendly femoral-groove design. TKA clinical-functional outcome and complications were highly influenced by the bearing geometry and kinematic pattern of prosthetic designs.

Evidence in focus

Study summary

Iriuchishima T, Ryu K. *J Knee Surg* (2017)*

Supporting healthcare
professionals for over 150 years

JOURNEY® II BCS restores function comparable to bicruciate retaining Oxford Unicompartmental Knees (UKA)

The bicruciate stabilising design reproduces anterior and posterior cruciate ligament function and native knee rollback

Study design

A single-surgeon retrospective study comparing outcomes at 6–9 months follow-up of patients treated with:

- JOURNEY II BCS: n=64 patients (mean age, 71.3 years $\pm$ 7.2 years)
- UKA: n=50 patients (mean age, 73.8 years $\pm$ 6 years)
- Control group: contralateral asymptomatic knees of subjects with UKA

Key results

Post-operative lateral knee radiograph (full flexion) showed:

- No significant difference in rollback ratio or knee flexion angle among the three groups (Table 1)
- Significant correlation between rollback ratio and knee flexion angle among the three groups ($p=0.002$) (Table 1)

	JOURNEY II BCS	UKA	Control
Rollback ratio %, mean (standard deviation)	37.9 ($\pm$ 4.9)	35.7 ($\pm$ 4.2)	35.3 ($\pm$ 4.8)
Flexion angle degrees, mean (standard deviation)	123.8 ($\pm$ 8.4)	125.4 ($\pm$ 7.5)	127 ($\pm$ 10.3)

Table 1. Rollback ratio and flexion angle measurement

Conclusion

- JOURNEY II BCS showed no significant difference in rollback ratio when compared with UKA or asymptomatic control knees.
- The implant design is likely to reproduce native anterior and posterior cruciate function and native knee rollback.

Considerations

- Pre-operative evaluation of rollback ratio and knee flexion was not performed

Study citation

*Iriuchishima T, Ryu K. A Comparison of Rollback Ratio between Bicruciate. Substituting Total Knee Arthroplasty and Oxford Unicompartmental Knee Arthroplasty. *J Knee Surg*. 2017 Jul 25. [Epub ahead of print]

64 JOURNEY™ II BCS vs 50 UKA vs 50 Healthy Knees

- No significant difference in rollback ratio or ROM among 3 groups
- Significant correlation between rollback ratio and ROM
- JOURNEY II BCS is likely to reproduce

Iriuchishima et al.

Etudes

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Nodzo, Mayman et al.

23° ROM improvement at 1-year follow-up vs standard PS design.¹

Kaneko et al.

Provides medio-lateral stability in the mid-flexion range and reproduces the **same medial pivot as normal knees**, aiding in patient recovery.²

Harris et al.

Significant improvements in patient-reported KSS and movement-associated pain at 24 months post-TKA.³

Mayman et al.

JOURNEY™ II BCS associated with **significantly reduced total hospital cost, less likely to have 30 days readmission and significantly reduced hospital stays** compared with other total knee arthroplasty (TKA) systems.⁴

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Hommel et al.

Significant improvements over baseline in Knee Score, Function Score and range of motion (ROM) with **patients averaging 127° ROM**.⁵

Iriuchishima et al.

Restores function comparable to bi-cruciate retaining Oxford Uni Knees (UKA).⁶

Grieco, Komistek et al.

Exhibits normal-like knee kinematic patterns suggesting the dual cam-post design and asymmetric articular geometries **adequately replicate ACL and PCL function**.⁷

Satisfaction patient majorée

Nodzo, Mayman et al.

Significantly better KSS scores than **standard PS** designs at one year follow-up, resulting in **reducing patient dissatisfaction**.¹

Kaneko et al.

Reduces medio-lateral instability in mid-flexion and **improve patient satisfaction**.²

Harris et al.

Increased patient satisfaction scores from 33 at 6-month follow-up to 35 at 24-month follow-up.³

¹ Nodzo SR, et al. *Tech Orthop*. 2018;33:37-41.

² Kaneko T, et al. *J Orthop*. 2017;14(1):201-206.

³ Harris AI, et al. *Arthroplasty Today*. 2018;4(2):240-243.

⁴ Mayman DJ, et al. Poster presented at: ISPOR Symposium; May 19-23, 2018; Ba

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